

# Dental




Elevators	02
Extracting Forceps	03
Tooth Forceps	06
Wire Cutting Forceps	08



## Bein Elevator

140mm long



### Straight


Code	Description
PKDEN/8348	2mm wide
PKDEN/8346	2.5mm wide
PKDEN/8344	3mm wide
PKDEN/8342	3.5mm wide
PKDEN/8340	4mm wide

### Curved

Code	Description
PKDEN/8352	2mm wide
PKDEN/8347	2.5mm wide
PKDEN/8345	3mm wide
PKDEN/8343	3.5mm wide
PKDEN/8341	4mm wide

## Cryer-White Elevator


150mm long



Code	Description
PKDEN/8261	Left
PKDEN/8263	Right

## Cryer Elevator

150mm long



Code	Description
PKPL/069	Right
PKPL/069a	Left


## Coleman Elevator

145mm long



Code	Description
PKDEN/8294	Left


Code	Description
PKDEN/8295	Right



Code	Description
PKDEN/8273	Straight

## Warwick-James Elevator


150mm long



Code	Description
PKDEN/8275	Straight
PKDEN/8249	Left
PKDEN/8251	Right

## Coupland Elevator



Curved, 150mm long



Code	Description
PKPL/068	3mm wide - Ash no1
PKPL/068a	4mm wide - Ash no2
PKPL/068b	5mm wide - Ash no3

## Winter Elevator

105mm long




Code	Description
PKDEN/8338	Left
PKDEN/8358	Right

# Extracting Forceps

## Extracting Forceps

For lower molars with decayed or broken crowns either side molars  
180mm long




Code	Description
PKDEN/8578	No 87




## Extracting Forceps

For lower molars either side  
165mm long




Code	Description
PKDEN/8522	No 22




## Extracting Forceps

For upper roots  
180mm long




Code	Description
PKDEN/8530	No 30




## Extracting Forceps

For upper roots and incisors  
165mm long



Code	Description
PKDEN/8522a	No 29




## Extracting Forceps

For upper molars  
180mm long




Code	Description
PKDEN/8588	No 89
PKDEN/8590	No 90




## Extracting Forceps

For upper laterals and canines  
170mm long



Code	Description
PKDEN/8500	No 1, 6mm wide
PKDEN/8502	No 1, 5mm wide




## Extracting Forceps

For lower molars  
165mm long




Code	Description
PKDEN/8520	No 21




## Extracting Forceps

For lower incisors and canines  
170mm long




Code	Description
SC-DEN/8504	No 4, 4mm wide



# Extracting Forceps

## Extracting Forceps


For lower premolars  
170mm long



Code	Description
PKDEN/8508	No 8, 4mm wide

## Extracting Forceps


For lower molars  
150mm long



Code	Description
PKDEN/8566	No 73, 3.5mm wide

## Extracting Forceps


For lower bicuspid and canines  
170mm long



Code	Description
PKDEN/8512	No 13, 4mm wide

## Extracting Forceps


For lower roots  
150mm long



Code	Description
SC-DEN/8568	No 74, 4mm wide

## Extracting Forceps


For upper premolar  
180mm long



Code	Description
PKDEN/8506	No 307, 5mm wide

## Extracting Forceps


For lower roots and premolars  
150mm long



Code	Description
PKDEN/8570	No 74n, 3mm wide

## Extracting Forceps


For upper molars  
180mm long



Code	Description
PKDEN/8516	No 17, right
PKDEN/8518	No 18, left

## Extracting Forceps

For lower bicuspid and canines  
150mm long




Code	Description
PKDEN/8572	No 75, 4mm wide

# Extracting Forceps

## Extracting Forceps

For lower molars with decayed or broken crowns either side  
150mm long




Code	Description
PKDEN/8576	No 86a




## Extracting Forceps

For upper roots, narrow  
185mm long




Code	Description
PKDEN/8552	No 51a, 3mm wide




## Extracting Forceps

For lower roots either side  
165mm long




Code	Description
PKDEN/8532	No 33, 3.5mm wide




## Extracting Forceps

For upper wisdoms either side  
185mm long



Code	Description
PKDEN/8560	No 67, 7mm wide




## Extracting Forceps

For upper roots  
185mm long



Code	Description
PKDEN/8550	No 51, 4mm wide




## Extracting Forceps

For lower wisdoms either side  
185mm long



Code	Description
PKDEN/8574	No 79, 5mm wide






## Cryer Tooth Forceps

For upper incisors and roots  
185mm long



Code
PKDEN/8010




## Cryer Tooth Forceps

For incisors and roots  
185mm long



Code	Description
PKDEN/8013	No 150a




## Cryer Tooth Forceps

For incisors and roots  
185mm long



Code	Description
PKDEN/8014	No 62




## Tooth Forceps

For lower incisors, bicuspid and roots  
185mm long



Code	Description
PKDEN/8036	No 203




## Cryer Tooth Forceps

For upper incisors and roots  
185mm long



Code	Description
PKDEN/8012	No 150




## Tooth Forceps

For lower molars either side  
185mm long




Code	Description
PKDEN/8039	No 23



# Tooth Forceps

## Tooth Forceps



For lower molars  
185mm long



Code	Description
PKDEN/8042	No 17

## Tooth Forceps


For lower roots, incisors and bicuspid  
175mm long



Code	Description
PKDEN/8045	No 103

## Tooth Forceps

For lower wisdom teeth  
185mm long



Code	Description
PKDEN/8043	No 222

## Tooth Forceps for Children

160mm long




Code	Description
PKDEN/8037	No 151as, for lower teeth

Code	Description
PKDEN/8017	No 150as, for upper teeth

## Tooth Forceps



For lower molars either side  
175mm long



Code	Description
PKDEN/8040	No 16

## Cryer Tooth Forceps for Children

For upper incisors, canines and bicuspid,  
160mm long




Code	Description
PKDEN/8011	No 1505

Code	Description
PKDEN/8033	No 1515

# Wire Cutting Forceps

## Top Cutting Forceps


Tungsten carbide edge, for bone pins up to 2.5mm  $\varnothing$   
gold handles



Code	Description
PKOR/06	215mm long

## Side Cutting Forceps

Tungsten carbide edge, for bone pins up to 2.2mm  $\varnothing$   
gold handles



Code	Description
PKOR/332	230mm long

## Top Cutting Forceps

Tungsten carbide edge, for arthrodesis wire up to 1.5mm  $\varnothing$   
gold handles



Code	Description
PKOR/06a	145mm long