




ENT

| | |
|--------------------------------------|----|
| Middle Ear Instruments | 02 |
| Picks and Scoops | 04 |
| Hooks and Probes | 05 |
| Retractors | 06 |
| Suction Tubes | 07 |
| Speculum | 08 |
| Aural Forceps | 10 |
| Artery Forceps | 12 |
| Tonsil, Turbinate and Septum Forceps | 13 |
| Dressing Forceps | 14 |
| Other Forceps | 15 |
| Mouth Gags | 16 |
| Tongue Plates and Depressors | 17 |
| Punches | 18 |
| Rongeurs | 22 |
| Scissors | 24 |
| Gouges, Knives and Elevators | 26 |
| Curettes and Myringotomes | 27 |
| Elevators and Dissectors | 28 |
| Miscellaneous | 29 |

Middle Ear Instruments

Ironside Instruments


160mm long, black finish



| Code | Description |
|------------|---|
| PVCENT200a | Knife 1, vertical incision |
| PVCENT200b | Knife 2, horizontal incision |
| PVCENT200c | Elevator |
| PVCENT200d | Needle point |
| PVCENT200e | Hook, superior elevation, angle 45° |
| PVCENT200f | Hook, inferior elevation, angle 45° |
| PVCENT200g | Hook, anterior elevation, angle 90° right |
| PVCENT200h | Hook, anterior elevation, angle 90° left |
| PVCENT200i | Measuring device |

Rosen Instruments


160mm long, black finish



| Code | Description |
|------------|--------------------------------------|
| PVCENT202a | Curette, tympanic ring, small |
| PVCENT202b | Curette, tympanic ring, large |
| PVCENT202c | Probe, footplate movement test, fine |
| PVCENT202d | Mobiliser, slight curve |
| PVCENT202e | Mobiliser, full curve |
| PVCENT202f | Scoop, size 000 |
| PVCENT202g | Scoop, size 00 |
| PVCENT202h | Scoop, size 0 |
| PVCENT202i | Scoop, size 1 |
| PVCENT202j | Fenestrator, oval window |

Schuknecht Instruments


160mm long, black finish



| Code | Description |
|------------|---------------------------|
| PVCENT201a | Hook upward, 0.4mm |
| PVCENT201b | Hook downward, 0.4mm |
| PVCENT201c | Hook upward, 0.7mm |
| PVCENT201d | Hook downward, 0.7mm |
| PVCENT201e | Elevator, oval |
| PVCENT201f | Hook, right-angled upward |
| PVCENT201g | Pick, straight |
| PVCENT201h | Spatula |

Beale Instruments


165mm long, satin finish



| Code | Description |
|------------|-------------------------------------|
| PVCENT204a | Tympanoplasty probe, fine blunt end |
| PVCENT204b | Raspatory, slight curve |
| PVCENT204c | Raspatory, slight curve, small |
| PVCENT204d | Rugine, double ended, sharp |
| PVCENT204e | Rugine, double ended |
| PVCENT204f | Spatula, graft smoothing |

Beale Instruments


170mm long, satin finish



| Code | Description |
|------------|--|
| PVCENT203a | Knife, Stapedius tendon cutting |
| PVCENT203b | Raspatory, mucous removal |
| PVCENT203c | Osteotome, footplate |
| PVCENT203d | Perforator, elevation fixed footplate |
| PVCENT203e | Hook, removal of footplate, angled downward |
| PVCENT203f | Hook, removal of footplate, angled upward |
| PVCENT203g | Hook, bone removal from oval window, angled left |
| PVCENT203h | Applicator, vein graft for right ear |
| PVCENT203i | Applicator, vein graft for left ear |
| PVCENT203j | Curette, removal bony rim |
| PVCENT203k | Curette |

Hough Instruments


165mm long, stainless steel



| Code | Description |
|------------|--|
| PVCENT205a | Excavator, "whirlybird", right, large |
| PVCENT205b | Excavator, "whirlybird", left, large |
| PVCENT205c | Excavator, "whirlybird", right, small |
| PVCENT205d | Excavator, "whirlybird", left, small |
| PVCENT205e | Pick, stapedectomy footplate, angled 90°, 0.33mm |
| PVCENT205f | Pick, stapedectomy footplate, angled 45°, 0.5mm |
| PVCENT205g | Pick, stapedectomy footplate, angled 15°, 0.5mm |
| PVCENT205h | Hook, "hoe", excavator, angled 90° |
| PVCENT205i | Hook, "hoe", excavator, angled 45° |
| PVCENT205j | Pick, curved |

Cawthorne Hooks

165mm long, black finish



| Code | Description |
|------------|------------------|
| PVCENT206a | Slight curve |
| PVCENT206b | Medium curve |
| PVCENT206c | Full curve |
| PVCENT206d | Double curve |
| PVCENT206e | 90° angle, 0.5mm |
| PVCENT206f | 90° angle, 1.0mm |
| PVCENT206g | 45° angle, 0.5mm |
| PVCENT206h | Retrograde |
| PVCENT206i | Straight |

Protect your delicate ENT instruments with our custom baskets



Picks and Scoops

Musgrove Cerumen Pick and Wool Carrier

Double ended



| Code | Description |
|-----------|-------------|
| PVCENT024 | 178mm long |

Formby Cerumen Scoop and Hook

Double ended



| Code | Description |
|-----------|-------------|
| PVCENT026 | 178mm long |

Hovell Loop


Fine, for removing hard wax



| Code | Description |
|-----------|-------------|
| PVCENT027 | 152mm long |

Lempert Scoop

Oval, 200mm long



| Code | Description |
|------------|-------------|
| PVCENT066a | 2mm x 2.5mm |
| PVCENT066b | 2mm x 3.5mm |
| PVCENT066c | 2mm x 4mm |
| PVCENT066d | 2mm x 4.5mm |
| PVCENT066e | 3mm x 5.5mm |
| PVCENT066f | 3mm x 6mm |

House Curette

Double ended, 180mm long



| Code | Description |
|------------|-------------|
| PVCSH5213 | Straight |
| PVCSH5213a | Angled |

Balance Scoop and Curette



| Code | Description |
|-----------|-------------|
| PVCENT038 | 230mm long |

Hartmann Hook

Sharp, straight shaft



| Code | Description |
|-----------|-------------|
| PVCENT021 | 140mm long |


Gardiner Brown Hook and Imray Loop

Double ended



| Code | Description |
|-----------|-------------|
| PVCENT025 | 127mm long |

Cawthorne Cerumen Hook



| Code | Description |
|-----------|-------------|
| PVCENT023 | 143mm long |

St Bart's Hook


Blunt, straight handle



| Code | Description |
|-----------|-------------|
| PVCENT022 | 178mm long |

Tracheal Hook

Sharp



| Code | Description |
|--------------|-------------|
| PVCENTSH5672 | 165mm long |

McIndoe Skin Hook


Fine, with flat handle



| Code | Description |
|-----------|-------------|
| PVCENT057 | 191mm long |

Tracheal Single Hook Retractor

Sharp



| Code | Description |
|-----------|-------------|
| PVCSH5671 | 152mm long |

Cawthorne Probe

Silver, flattened end



| Code | Description |
|-----------|-------------|
| PVCENT019 | 140mm long |

Tracheal Double Hook

Blunt



| Code | Description |
|--------------|-------------|
| PVCENTSH5673 | 105mm long |

Jobson Horne Probe

Double ended, serrated tip and smooth ring




| Code | Description |
|------------|-------------|
| PVCENT020a | 140mm long |
| PVCENT020b | 178mm long |

Retractors


Mollison Retractor

Self retaining, curved



| Code | Description |
|------------|------------------------------|
| PVCENT034a | 4-4 teeth, adult, 165mm long |
| PVCENT034b | 2-2 teeth, child, 133mm long |


Wallace Modified Plester Retractor



| Code | Description |
|-----------|-------------------------|
| PVCENT062 | 2 x 2 teeth, 140mm long |

Plester-Jansen Retractor


Solid blade one side,
140mm long



| Code | Description |
|------------|-------------|
| PVCENT061A | Left |
| PVCENT061B | Right |
| PVCENT061C | 3 x 3 Teeth |

Shah Intra Aural Retractor

Satin finish
152mm long



| Code | Description |
|------------|------------------------------|
| PVCENT060a | 4-4 teeth, adult, 165mm long |
| PVCENT060b | Third blade |
| PVCENT060c | Intra aural blade |
| PVCENT060d | Post aural blade |

Hurd Pillar Retractor and Tonsil Dissector



| Code | Description |
|-----------|-------------|
| PVCENT147 | 197mm long |

Mollison Anterior Pillar Retractor and Tonsil Enucleator

Semi sharp



| Code | Description |
|-----------|-------------|
| PVCENT148 | 210mm long |

Ipswich Pattern Sucker



| Code | Description |
|------------|-------------|
| PVCENT050a | 267mm long |

Lempert Suction and Irrigation Tubes


With tucker valve and stilette for cleaning stainless steel, 203mm



| Code | Description |
|------------|-------------|
| PVCENT044 | 6F.G. 2mm |
| PVCENT044a | 9F.G. 3mm |

Formby Suction Tube

Angled, conical, tip size 10 F.G. stainless steel



| Code | Description |
|-----------|-------------|
| PVCENT042 | 210mm long |

Magill Suction Tubes


Set of four with conical ends, one curved, three angled, handle, stainless steel, 210mm



| Code | Description |
|------------|---------------|
| PVCENT045 | Various sizes |
| PVCENT046a | 3.0mm ø |
| PVCENT046b | 3.5mm ø |
| PVCENT046c | 4.0mm ø |

Barron Suction Tube

Small size, 4 F.G. Not for use with Zoellner ends



| Code | Description |
|-----------|-------------|
| PVCENT041 | 229mm long |

Zoellner Suction Tube

With stilette, stainless steel




| Code | Description |
|-----------|-------------------|
| PVCENT047 | 170mm long, 2mm ø |




House Suction and Irrigation Tubes

Stainless steel, 152mm long



| Code | Description |
|------------|------------------|
| PVCENT043 | 8 F.G. x 6 F.G. |
| PVCENT043a | 10 F.G. x 7 F.G. |
| PVCENT043b | 12 F.G. x 8 F.G. |

Rowbotham Suction Tube



| Code | Description |
|-----------|-------------|
| PVCENT169 | 286mm long |

Speculum

Rosen Slotted Speculum

Set of four, round ends, black finish
38mm long


| Code | Description |
|------------|-------------|
| PVCENT049j | Set of 4 |
| PVCENT049f | 4mm |
| PVCENT049g | 5mm |
| PVCENT049h | 6mm |
| PVCENT049i | 7mm |



Tumarkin Speculum

Set of eight, tapered slot, black finish
50mm long

| Code | Description |
|------------|-------------|
| PVCENT052 | Set of 8 |
| PVCENT053a | 4.5mm ø |
| PVCENT053b | 6.5mm ø |
| PVCENT053c | 7.0mm ø |
| PVCENT053d | 7.5mm ø |
| PVCENT053e | 8.0mm ø |
| PVCENT053f | 8.5mm ø |
| PVCENT053g | 9.0mm ø |
| PVCENT053h | 11.0mm ø |



Rosen Aural Speculum

Round, 40mm long

| Code | Description |
|------------|-------------|
| PVCENT049a | 3mm |
| PVCENT049b | 4mm |
| PVCENT049c | 5mm |
| PVCENT049d | 6mm |
| PVCENT049e | 7mm |



Heath Speculum

Set of four, black finish
50mm long

| Code | Description |
|------------|-------------|
| PVCENT207e | Set of 4 |
| PVCENT207a | 4.0mm ø |
| PVCENT207b | 5.0mm ø |
| PVCENT207c | 6.0mm ø |
| PVCENT207d | 7.0mm ø |




Shea Speculum

Set of five, oblique ends, black finish, 38mm long

| Code | Description |
|------------|-------------|
| PVCENT050 | Set of 5 |
| PVCENT051 | 3mm |
| PVCENT051a | 4mm |
| PVCENT051b | 5mm |
| PVCENT051c | 6mm |
| PVCENT051d | 7mm |




Peters Seigle Speculum



| Code | Description |
|------------|--|
| PVCENT366 | Peters ear speculum black chrome finish set of 5, sizes 0 to 4 |
| PVCENT189i | 2.2mm x 3.2mm |
| PVCENT189j | 2.8mm x 4.2mm |
| PVCENT189k | 4.3mm x 5.8mm |
| PVCENT189l | 5.2mm x 7.6mm |
| PVCENT189m | 6.0mm x 7.6mm |
| PVCENT189d | 6mm |
| PVCENT189e | 7mm |
| PVCENT189f | 8mm |
| PVCENT189g | 9mm |
| PVCENT189h | 10mm |

St Clair Thompson Nasal Speculum


| Code | Description |
|------------|-------------|
| PVCENT072 | 38mm long |
| PVCENT072a | 51mm long |
| PVCENT072b | 64mm long |
| PVCENT072c | 89mm long |



Killian Nasal Speculum

With Stop

| Code | Description |
|------------|-------------|
| PVCENT081a | 51mm blade |
| PVCENT081b | 64mm blade |
| PVCENT081c | 76mm blade |
| PVCENT081d | 89mm blade |




Without Stop

| Code | Description |
|-------------|-------------|
| SC-ENT/082b | 76mm blade |

Thudichum Nasal Speculum

Set of seven, sizes 1-7,
70mm long

| Code | Description |
|------------|-------------|
| PVCENT073a | Standard |
| PVCENT073b | Insulated |



Standard

| Code | Description |
|------------|---------------------|
| PVCENT073a | Size 1 – 12mm blade |
| PVCENT073a | Size 2 – 14mm blade |
| PVCENT073a | Size 3 – 16mm blade |
| PVCENT073a | Size 4 – 18mm blade |
| PVCENT073a | Size 5 – 20mm blade |
| PVCENT073a | Size 6 – 22mm blade |
| PVCENT073a | Size 7 – 24mm blade |
| PVCENT073a | Size 8 – 27mm blade |

Insulated


| Code | Description |
|------------|---------------------|
| PVCENT073b | Size 1 – 12mm blade |
| PVCENT073b | Size 2 – 14mm blade |
| PVCENT073b | Size 3 – 16mm blade |
| PVCENT073b | Size 4 – 18mm blade |
| PVCENT073b | Size 5 – 20mm blade |
| PVCENT073b | Size 6 – 22mm blade |
| PVCENT073b | Size 7 – 24mm blade |
| PVCENT073b | Size 8 – 27mm blade |

Lempert Endaural Speculum



| Code | Description |
|------------|-------------|
| PVCENT060e | 133mm long |


Pilcher Nasal Speculum



| Code | Description |
|-----------|-------------|
| PVCENT071 | 152mm |


Cottle Nasal Speculum

127mm long



| Code | Description |
|------------|-------------|
| PVCENT192c | 20mm blade |
| PVCENT192b | 32mm blade |
| PVCENT192 | 41mm blade |
| PVCENT192a | 60mm blade |

Consistent Quality Inspection Process



Get in touch to query your instrument requirements at vanite@cwccsa.com


Aural Forceps



Triangular Jaw



Smooth Jaw



Serrated Jaw



Tenaculum Jaw



Cup Jaw



Grunwald Punch



Henckel Jaw

Beale Aural Forceps

Micro, black finish, triangular flat jaws, serrated



| Code | Description |
|-----------|------------------|
| PVCENT183 | 70mm to shoulder |

Hough Aural Forceps

Micro, extra fine oval cup jaws, black finish



| Code | Description |
|-----------|------------------|
| PVCENT059 | 70mm to shoulder |

Cawthorne Aural Forceps

Micro, black finish, extra fine
70mm to shoulder



| Code | Description |
|-----------|---------------|
| PVCENT182 | Smooth jaws |
| PVCENT199 | Serrated jaws |

Ormerod Aural Forceps


Micro, extra fine, black finish
70mm to shoulder



| Code | Description |
|------------|------------------------------|
| PVCENT184a | Serrated, angled upwards |
| PVCENT184b | 1mm cup jaw |
| PVCENT184c | Elongated cup jaw, 0.5 x 4mm |

Hartmann Aural Forceps

Micro, black finish 5mm to shoulder



| Code | Description |
|------------|-------------------------------|
| PVCENT058 | Serrated jaws |
| PVCENT058a | Tenaculum jaws |
| PVCENT058b | Cup shape jaws |
| PVCENT058c | Curved left, serrated |
| PVCENT058d | Curved right, serrated |
| PVCENT058e | Angled up, serrated |
| PVCENT058f | Serrated jaws, satin finish |
| PVCENT058s | Pre-sterile - single use only |

Shea Aural Forceps


Micro, oval cup jaw, black finish
70mm to shoulder



| Code | Description |
|------------|-----------------|
| PVCENT069 | Standard |
| PVCENT069a | Curved to left |
| PVCENT069b | Curved to right |
| PVCENT069c | Curved upward |

Hartmann Aural Forceps


Standard, satin finish
75mm to shoulder



| Code | Description |
|------------|-----------------------|
| PVCENT181 | Serrated jaws |
| PVCENT181a | Tenaculum jaws |
| PVCENT181b | Cup shape jaws |
| PVCENT181c | Grunwald punch action |

Heath Aural Forceps


Light, with lifting discs, serrated jaws



| Code | Description |
|-----------|------------------|
| PVCENT175 | 83mm to shoulder |

Henckel Aural Forceps


Fine fenestrated scoop blades



| Code | Description |
|-----------|------------------|
| PVCENT176 | 75mm to shoulder |

Struycken Cutting Forceps

75mm to shoulder




| Code | Description |
|-----------|-----------------------|
| PVCENT178 | Fenestrated under jaw |

Artery Forceps

Birkett Tonsil Artery Forceps


Fine, 178mm long



| Code | Description |
|------------|-------------|
| PVCENT151a | Straight |
| PVCENT151b | Curved |


Negus Tonsil Artery Forceps

180mm long



| Code | Description |
|------------|-------------|
| PVCENT197 | Large curve |
| PVCENT197a | Small curve |


Denis Browne Tonsil Artery Forceps



| Code | Description |
|-----------|-------------|
| PVCENT155 | 216mm long |

Sawtell Tonsil Artery Forceps


178mm long



| Code | Description |
|-----------|-------------|
| PVCENT152 | Large curve |
| PVCENT153 | Small curve |

Fickling Artery Forceps


190mm long



| Code | Description |
|-----------|-------------|
| PVCPLO44b | 1-2 teeth |
| PVCPLO44a | Serrated |

Scott Tonsil Artery Forceps

Curved to side, screw joint




| Code | Description |
|-----------|-------------|
| PVCENT154 | 180mm long |

Tonsil, Turbinate and Septum Forceps

Denis Browne Tonsil Holding Forceps

Angled shanks, no rack



| Code | Description |
|-----------|-------------|
| PVCENT144 | 200mm long |

Jansen-Middleton Nasal Septum Forceps


Compound action
178mm long



| Code | Description |
|----------|-------------|
| PVCPLO32 | 5mm jaw |


Irwin Moore Nasal Turbinate Forceps

191mm long




| Code | Description |
|------------|-------------|
| PVCENT077d | 5mm jaw |

Walsham Septum Forceps



| Code | Description |
|-----------|----------------------|
| PVCPLO30a | Straight, 215mm long |
| PVCPLO30b | Right, 215mm long |
| PVCPLO30c | Left, 215mm long |

Luc Nasal Turbinate Forceps



Heart shaped, 178mm long


| Code | Description |
|------------|----------------|
| PVCENT078a | Small 6mm ø |
| PVCENT078b | Medium 7.5mm ø |
| PVCENT078c | Large 9.0mm ø |

Oval, 191mm long

| Code | Description |
|------------|----------------|
| PVCENT078d | Small 6mm ø |
| PVCENT078e | Medium 7.5mm ø |
| PVCENT078f | Large 9.0mm ø |

Hartmann Conchotomes

Circular cutting 102mm to shoulder




| Code | Description |
|------------|-------------|
| PVCENT080a | 5mm ø |
| PVCENT080b | 7mm ø |
| PVCENT080c | 9mm ø |
| PVCENT080d | 11mm ø |

Dressing Forceps

Hartmann Dressing Forceps


Serrated jaw, screw joint



| Code | Description |
|-----------|-------------------|
| PVCENT076 | 108mm to shoulder |

Wilde Dressing Forceps

Angular, serrated jaws



| Code | Description |
|-----------|-------------------|
| PVCENT075 | 102mm to shoulder |

Tilley Dressing Forceps


Straight, fine, serrated ends



| Code | Description |
|-----------|----------------------------|
| PVCENT077 | 2mm wide, 76mm to shoulder |

Irwin Moore Thyroid Clamp


Thin jaws, angled to side



| Code | Description |
|-----------|-------------|
| PVCENT170 | 165mm long |




Trousseau Tracheal Dilator



| Code | Description |
|------------|------------------|
| PVCENT172 | Adult 140mm long |
| PVCENT172a | Child 125mm long |

Wilson Bone Nibbling Forceps


Fine, curved, satin finish
152mm long



| Code | Description |
|-----------|-----------------------------|
| PVCENT054 | Downward cutting - 3mm jaws |
| PVCENT055 | Upward cutting - 3mm jaws |


Cawthorne Nibbling Forceps

Narrow, 3mm jaws, wide handle



| Code | Description |
|-----------|-------------|
| PVCENT068 | 159mm long |

Denis Browne Collar Crushing Forceps



| Code | Description |
|------------|-------------|
| PK01SH4062 | 140mm long |

Lempert Nibbling Forceps


152mm long



| Code | Description |
|------------|----------------------|
| PVCENT067a | Straight, 2.5mm jaws |
| PVCENT067b | Angled, 3.0mm jaws |

Mouth Gags


Doyen Mouth Gag



| Code | Description |
|-----------|-------------|
| PVCENT88d | Adult |
| PVCENT88e | Child |

Boyle-Davis Mouth Gag


Frame only



| Code | Description |
|------------|-------------|
| PVCENT083a | Adult |
| PVCENT083b | Child |


Kilner Dott Mouth Gag

Frame only




| Code | Description |
|------------|---|
| PVCENT087a | Adult |
| PVCENT088e | Set, including frame, small, medium, and large blades |

Fergusson Ackland Mouth Gag




| Code | Description |
|------------|--------------|
| PVCENT209 | Adult, 210mm |
| PVCENT209a | Child, 155mm |

Mason Mouth Gag



| Code | Description |
|-----------|--------------|
| PVCENT210 | Adult, 210mm |

Draffin Bipod



Adult

| Code | Description |
|------------|---|
| PVCENT083 | Adult, 4 hole x 480mm long half circle end for use with a drape |
| PVCENT083E | Adult, 4 hole x 480mm long with round ends to be used with Magauran plate |
| PVCENT083E | Adult, 4 hole x 480mm long with round ends to be used with Magauran plate |

Child


| Code | Description |
|------------|---|
| PVCENT083C | Child, 3 hole x 425mm long half circle end for use with a drape |
| PVCENT083D | Child, 3 hole x 425mm long with round ends to be used with Magauran plate |
| PVCENT083F | Child, 4 hole x 425mm long with round ends to be used with Magauran plate Alder Hey Children's Hospital pattern |
| PVCENT083D | Child, 3 hole x 425mm long with round ends to be used with Magauran plate |
| PVCENT083F | Child, 4 hole x 425mm long with round ends to be used with Magauran plate Alder Hey Children's Hospital pattern |
| PK01SH5104 | Magauran plate 480mm - for use with Draffin bipods |

Magauran Plate

| Code | Description |
|------------|-------------------------------------|
| PK01SH5104 | 480mm - for use with Draffin bipods |

Tongue Plates and Depressors

Doughty Tongue Plate



Plain

| Code | Description |
|------------|-------------------------|
| PVCENT084a | 65mm - 1.6mm-1.9mm slot |
| PVCENT084b | 75mm - 1.6mm-1.9mm slot |
| PVCENT084c | 90mm - 4mm-4.5mm slot |
| PVCENT084d | 100mm - 4mm-4.5mm slot |
| PVCENT084e | 115mm - 4mm-4.5mm slot |

Slotted

| Code | Description |
|------------|-------------------------|
| PVCENT085a | 65mm - 1.6mm-1.9mm slot |
| PVCENT085b | 75mm - 1.6mm-1.9mm slot |
| PVCENT085c | 90mm - 4mm-4.5mm slot |
| PVCENT085d | 100mm - 4mm-4.5mm slot |
| PVCENT085e | 115mm - 4mm-4.5mm slot |

With anaesthetic tube

| Code | Description |
|------------|-------------------------|
| PVCENT086a | 65mm - 1.6mm-1.9mm slot |
| PVCENT086b | 75mm - 1.6mm-1.9mm slot |
| PVCENT086c | 90mm - 4mm-4.5mm slot |
| PVCENT086d | 100mm - 4mm-4.5mm slot |
| PVCENT086e | 115mm - 4mm-4.5mm slot |


Tongue Depressor

Flat stainless steel



| Code | Description |
|-----------|-------------|
| PVCENT143 | 184mm long |

Lack Tongue Depressor



| Code | Description |
|------------|-----------------------|
| PVCENT142a | Child, 13mm wide |
| PVCENT142b | Adolescent, 19mm wide |
| PVCENT142c | Adult, 22mm wide |

Kilner Dott Tongue Plate


With endotracheal tube groove



| Code | Description |
|------------|-------------|
| PVCENT088a | Small |
| PVCENT088b | Medium |
| PVCENT088c | Large |
| PVCENT088d | Extra small |

Punches


Nasal Cutting Forceps



Myles

Through cutting, 120mm to shoulder


| Code | Description |
|--------------|------------------|
| PVCENTSH5889 | 12.5 x 3.5mm jaw |



Gruenwald-Henke

Through cutting, 120mm to shoulder


| Code | Description |
|--------------|---------------|
| PVCENTSH5281 | 7 x 2mm jaw |
| PVCENTSH5282 | 9 x 2.5mm jaw |



Struycken


Through cutting, 120mm to shoulder

| Code | Description |
|--------------|--------------|
| PVCENTSH5280 | 11 x 1mm jaw |



Gruenwald Nasal Punch Forceps


108mm to shoulder



| Code | Description |
|------------|--|
| PVCENT077a | A fenestrated jaw, 3.5mm |
| PVCENT077b | Fenestrated lower jaw, narrow, 2.5mm |
| PVCENT077c | Fenestrated and tapered lower jaw, 2-3mm |

Seymour Jones Antrum Punch Forceps


159mm long



| Code | Description |
|-----------|-------------|
| PVCENT132 | 4mm jaw |

Backbiter


180mm long



| Code | Description |
|-------------|-------------|
| PK01SH4765b | 7 x 3mm jaw |

Backbiter


200mm long



| Code | Description |
|------------|-------------|
| PVCENT5392 | 7 x 2mm jaw |

Stamberger Rhinoforcep Antrium Punch


7mm x 3mm, 130mm long



| Code | Description |
|------------|----------------|
| PVCENT234 | Left back cut |
| PVCENT234a | Right back cut |

Tilley Henckel Ethmoid Punch


Graduated shank, 114mm to shoulder



| Code | Description |
|------------|-------------|
| PVCENT138a | 4.5mm |
| PVCENT138b | 5.5mm |
| PVCENT138c | 6.5mm |
| PVCENT138d | 8mm |

Ostrom Antrum Punch


Forcep, 178mm to shoulder



| Code | Description |
|-----------|-------------|
| PVCENT135 | 3mm jaw |

Mackay Modified Stamberger Punch

200mm long



| Code | Description |
|------------|-------------|
| PVCENT5392 | 7 x 2mm jaw |

Paterson Biopsy Forceps

Elongated cup, 3mm jaw




| Code | Description |
|------------|-------------|
| PVCENT174a | 305mm long |
| PVCENT174b | 355mm long |
| PVCENT174c | 406mm long |
| PVCENT174d | 508mm long |

Punches

Chevalier Jackson Endoscope Forcep

With crocodile action and serrated jaw


| Code | Description |
|-----------|------------------------------|
| PVCENT241 | Effective shaft length 300mm |
| PVCENT242 | Effective shaft length 400mm |
| PVCENT243 | Effective shaft length 500mm |
| PVCENT244 | Effective shaft length 600mm |



Chevalier Jackson Vocal Nodule Endoscope Forcep

With crocodile action and curved cup

| Code | Description |
|-----------|------------------------------|
| PVCENT245 | Effective shaft length 300mm |
| PVCENT246 | Effective shaft length 400mm |
| PVCENT247 | Effective shaft length 500mm |
| PVCENT248 | Effective shaft length 600mm |



Chevalier Jackson Vocal Nodule Endoscope Forcep

With crocodile action and straight cup


| Code | Description |
|-----------|------------------------------|
| PVCENT249 | Effective shaft length 300mm |
| PVCENT250 | Effective shaft length 400mm |
| PVCENT251 | Effective shaft length 500mm |
| PVCENT252 | Effective shaft length 600mm |



St Bartholomew's Ethmoid Punch


Scoop blades, 3mm bite, 114mm to shoulder

| Code | Description |
|------------|------------------|
| PVCENT140a | Straight |
| PVCENT140b | Upward cutting |
| PVCENT140c | Downward cutting |



Ferris Smith Kerrison Punch


Heavy, 5mm bite,
229mm long



| Code | Description |
|------------|------------------|
| PVCENT137a | Upward cutting |
| PVCENT137B | Downward cutting |

Citelli Antrum Punch


76mm to shoulder



| Code | Description |
|------------|-------------|
| PVCENT133a | 3mm jaw |
| PVCENT133b | 5mm jaw |

Hajek Antrum Punch


197mm long



| Code | Description |
|-----------|-------------|
| PVCENT134 | 5mm jaw |

Hajek Sphenoidal Punch

146mm to shoulder



Upward cutting

| Code | Description |
|------------|-------------|
| PVCENT141a | 4mm bite |
| PVCENT141b | 5mm bite |

Downward cutting

| Code | Description |
|------------|-------------|
| PVCENT141c | 4mm bite |
| PVCENT141d | 5mm bite |

Rongeurs

Beck Rongeur

Straight, for intervertebral disc



| Code | Description |
|-----------|-----------------------|
| PVCNEU125 | 2mm bite, 130mm shaft |

Pennybacker Rongeur

120mm shaft




| Code | Description |
|-----------|-------------|
| PVCNEU133 | 3mm bite |
| PVCNEU107 | 4mm bite |

Cushing Rongeur

130mm shaft



| Code | Description |
|-------------|-------------|
| PVCNEU127-2 | 2mm bite |
| PVCNEU127-3 | 3mm bite |
| PVCNEU127-4 | 4mm bite |
| PVCNEU127-5 | 5mm bite |
| PVCNEU127-6 | 6mm bite |



St Bart's Rongeur



114mm shaft



| Code | Description |
|------------|-------------|
| PVCENT140a | 3mm bite |

Cushing Rongeur

130mm shaft



150° upward cutting


| Code | Description |
|------------|-------------|
| PVCNEU128c | 2mm bite |
| PVCNEU128 | 3mm bite |
| PVCNEU128a | 4mm bite |
| PVCNEU128b | 5mm bite |

150° downward cutting

| Code | Description |
|------------|-------------|
| PVCNEU130c | 2mm bite |
| PVCNEU130 | 3mm bite |
| PVCNEU130a | 4mm bite |
| PVCNEU130b | 5mm bite |

St Bart's Rongeur

140°, 3mm, 114mm shaft



| Code | Description |
|-------------|------------------|
| PVCENT140b | Upward cutting |
| PVCENT140ac | Downward cutting |

Mackay Rongeur


140°, straight, 114mm shaft



| Code | Description |
|------------|-------------|
| PK01SH4249 | 5.5mm |

Weil-Blakesley Rongeur

110mm shaft



140° upward cutting


| Code | Description |
|------------|-------------|
| PVCENT194b | 3.6mm bite |
| PVCENT194e | 4.2mm bite |
| PVCENT194f | 4.8mm bite |

140° downward cutting

| Code | Description |
|------------|-------------|
| PVCENT194g | 3.6mm bite |
| PVCENT194h | 4.2mm bite |
| PVCENT194i | 4.8mm bite |

Mackay Rongeur


With suction tube



| Code | Description |
|-------------|-------------------------|
| PK01SH4249a | 3.5mm bite, 114mm shaft |

Weil-Blakesley Rongeur


120mm shaft



| Code | Description |
|------------|-------------|
| PVCENT193b | 3.0mm bite |
| PVCENT193 | 3.6mm bite |
| PVCENT194 | 4.2mm bite |
| PVCENT194a | 4.8mm bite |
| PVCENT194c | 5.8mm bite |
| PVCENT194d | 6.6mm bite |

Mackay Rongeur

With suction tube, 140° upward cutting



| Code | Description |
|-------------|-------------------------|
| PK01SH4249b | 3.5mm bite, 114mm shaft |

Scissors

Beckmann Turbinectomy Scissors


Curved



| Code | Description |
|-----------|-------------|
| PVCENT097 | 165mm long |


Fomon Scissors

Straight, serrated blades, angled handles



| Code | Description |
|-----------|-------------|
| PVCENT200 | 145mm long |

Scott Stevenson Tonsil Scissors



| Code | Description |
|-----------|-------------|
| PVCENT159 | 216mm long |

Cottle Scissors

Straight blades, angled handles, blunt/blunt



| Code | Description |
|-------------|-------------|
| PVCENTOK201 | 165mm long |

Wood Tonsil Scissors


Sharpened outer edges



| Code | Description |
|-----------|-------------|
| PVCENT156 | 178mm long |

Cottle Bulldog Scissors


Blunt/blunt, 115mm long



| Code | Description |
|-----------|-------------|
| PVCENT202 | 115mm long |

Wheeler Tonsil Scissors


Triangular form blades, serrated shanks



| Code | Description |
|-----------|-------------|
| PVCENT157 | 184mm long |

Joseph Scissors


Sharp/sharp, 150mm long



| Code | Description |
|-----------|-------------|
| PVCENT203 | Straight |
| PVCENT204 | Curved |

Tilley Tonsil Scissors


Curved



| Code | Description |
|-----------|-------------|
| PVCENT158 | 190mm long |


Heymann Turbinectomy Scissors

Angular



| Code | Description |
|-----------|-------------|
| PVCENT098 | 197mm long |


Wilson Tonsil Scissors



| Code | Description |
|-----------|-------------|
| PVCENT160 | 203mm long |


Cawthorne Scissors

Micro, extra fine, 70mm to shoulder



| Code | Description |
|------------|-------------|
| PVCENT186c | Left |
| PVCENT186d | Right |
| PVCENT186e | Straight |

Thyroid Scissors



| Code | Description |
|-----------|-------------|
| PVCENT171 | 152mm long |

Ormrod Aural Scissors


Micro, extra fine, curved upwards



| Code | Description |
|-----------|------------------|
| PVCENT185 | 70mm to shoulder |

Knapp Scissors


Slight curve, blunt points, 114mm long




| Code | Description |
|------------|--------------|
| PVCENT187a | Satin finish |
| PVCENT187b | Black finish |

Portmann Aural Scissors

Micro, satin finish 70mm to shoulder




| Code | Description |
|------------|-----------------|
| PVCENT186a | Curved to right |
| PVCENT186b | Curved to left |



Gouges, Knives and Elevators

Ballenger Nasal Cartilage Knife

Revolving blade



Straight, 203mm long

| Code | Description |
|------------|-------------|
| PVCENT096a | 6mm wide |
| PVCENT096b | 8mm wide |

Cranked, 197mm long

| Code | Description |
|------------|-------------|
| PVCENT096c | 6mm wide |

Sickle Knife



| Code | Description |
|------------|-------------|
| PKBUPNS002 | 188mm long |
| PVCENT271 | 160mm long |

Jenkins Gouge

140mm long



| Code | Description | Code | Description |
|-------------|-------------|------------|-------------|
| SC-ENT/037 | 1mm | PVCENT037g | 8mm |
| SC-ENT/037a | 2mm | PVCENT037h | 9mm |
| SC-ENT/037b | 3mm | PVCENT037i | 10mm |
| SC-ENT/037c | 4mm | PVCENT037j | 11mm |
| SC-ENT/037d | 5mm | PVCENT037k | 12mm |
| SC-ENT/037e | 6mm | PVCENT037l | 13mm |
| SC-ENT/037f | 7mm | PVCENT037m | 14mm |

Killian Nasal Gouge

127mm long



| Code | Description |
|------------|-------------|
| SC-ENT/093 | 4mm |

Guy's Hospital Nasal Gouge

Shallow groove and flattened back, 127mm long



| Code | Description |
|-----------|-------------|
| PVCENT094 | 4mm |

Tilley Nasal Gouge

V shaped, 127mm long



| Code | Description |
|------------|-------------|
| SC-ENT/095 | 4mm |

Cottle Elevator


Double ended, black chrome finish
200mm long



| Code |
|-----------|
| PVCENT224 |

Hill Elevator

Double ended, curved to side, rugine end



| Code | Description |
|-----------|----------------------|
| PVCENT091 | 5mm wide, 222mm long |

Tumarkin Atticotomy Elevator

Double ended




| Code | Description |
|------------|-------------|
| SC-ENT/036 | 191mm long |

Curettes and Myringotomes

Buck Ear Curette

Angled, blunt, 146mm long



| Code | Description |
|-------------|-------------|
| SC-ENT/028 | 2mm |
| SC-ENT/028a | 2.5mm |
| SC-ENT/028b | 3mm |
| SC-ENT/028c | 3.5mm |
| SC-ENT/028d | 4mm |

Nasal Curette

Straight, 185mm long



| Code | Description |
|--------------|-------------|
| PVCENT628702 | 5 x 2mm |

Beckmann Adenoid Curette


216mm long



| Code | Description |
|------------|-------------|
| PVCENT167a | 8mm blade |
| PVCENT167b | 10mm blade |
| PVCENT167c | 12mm blade |
| PVCENT167d | 14mm blade |
| PVCENT167e | 16mm blade |
| PVCENT167f | 18mm blade |

Agnew (Politzer) Myringotome

190mm long



| Code | Description |
|-----------|--------------------|
| PVCENT031 | Shaft length 100mm |

Trautmann Myringotome


178mm long



| Code | Description |
|-----------|------------------|
| PVCENT032 | Upward cutting |
| PVCENT033 | Downward cutting |

St Clair Thomson Adenoid Curette

With cage, 235mm long




| Code | Description |
|------------|-------------|
| PVCENT168a | 8mm |
| PVCENT168b | 10mm |
| PVCENT168c | 12mm |
| PVCENT168d | 14mm |
| PVCENT168e | 16mm |
| PVCENT168f | 18mm |

Elevators and Dissectors

Freer Instruments

Single ended, 190mm long



Knife

| Code | Description |
|------------|-------------------------------|
| PVCENT065a | Curved to side |
| PVCENT065c | Rounded edge, 8mm blade' |
| PVCENT065d | Flat rounded blade |
| PVCENT065e | Rounded blade, curved on flat |

Elevator

| Code | Description |
|------------|-------------|
| PVCENT065f | Dull edged |
| PVCENT065h | Sharp edged |

Chisel

| Code |
|------------|
| PVCENT0651 |

St Clair Thompson Dissector

159mm



| Code | Description |
|------------|-------------|
| PVCENT089a | Sharp |
| PVCENT089b | Blunt |

Negus Aspiring Dissector


With 6mm olive, 210mm long



| Code | Description |
|------------|-------------|
| PVCENT090a | Insulated |
| PVCENT090b | Insulated |

Freer Separator and Elevator

Double ended



| Code | Description |
|----------------|---|
| PVCENT065b | Blunt/sharp, 180mm long |
| PVCENT065g | Sharp, knife curved to side, 180mm long |
| PVCENT065j | Blunt/sharp, 200mm long |
| PVCENT065j-sal | Blunt/sharp, 188mm long |

Syme Tonsil Dissector

Double ended, one end blunt, one semi sharp



| Code | Description |
|-----------|-------------|
| PVCENT149 | 170mm long |

Gwynne-Evans Tonsil Dissector


Double ended, with serrated end and blunt end 197mm long



| Code | Description |
|------------|-------------|
| PVCENT146a | 6mm wide |
| PVCENT146b | 8mm wide |

Woakes Eustachian Catheters

Silver, Luer, 140mm




| Code | Description |
|------------|-------------|
| PVCENT015 | 7 F.G. |
| PVCENT015a | 8 F.G. |
| PVCENT015b | 9 F.G. |

| Code | Description |
|------------|-------------|
| PVCENT015c | 10 F.G. |
| PVCENT015d | 11F.G. |

Kramer Eustachian Catheters

Silver, Luer, 140mm



| Code | Description |
|-------------|-------------|
| SC-ENT/014 | 7 F.G. |
| SC-ENT/014a | 8 F.G. |
| SC-ENT/014b | 9 F.G. |

| Code | Description |
|-------------|-------------|
| SC-ENT/014c | 10 F.G. |
| SC-ENT/014d | 11 F.G. |

Myles Naso-Antral Perforator

171mm long



| Code | Description |
|-------------|-------------|
| SC-ENT/103a | 16 F.G. |
| SC-ENT/103b | 19 F.G. |

Tilley Antrum Burrs


171mm long



| Code | Description |
|-------------|-------------|
| SC-ENT/131a | 6mm ø |
| SC-ENT/131b | 7mm ø |
| SC-ENT/131c | 8mm ø |

Laryngeal Mirror

Mirror only, 130mm long




| Code | Description |
|--------------|-------------|
| PVCENT191-6 | 6mm ø |
| PVCENT191-8 | 8mm ø |
| PVCENT191-10 | 10mm ø |
| PVCENT191-12 | 12mm ø |
| PVCENT191-14 | 14mm ø |
| PVCENT191-16 | 16mm ø |
| PVCENT191-18 | 18mm ø |

| Code | Description |
|--------------|-------------|
| PVCENT191-20 | 20mm ø |
| PVCENT191-22 | 22mm ø |
| PVCENT191-24 | 24mm ø |
| PVCENT191-26 | 26mm ø |
| PVCENT191-28 | 28mm ø |
| PVCENT191-30 | 30mm ø |

Laryngeal Mirror

Complete with metal handle, 260mm long



| Code | Description |
|--------------|-------------|
| PVCENT192-6 | 6mm ø |
| PVCENT192-8 | 8mm ø |
| PVCENT192-10 | 10mm ø |
| PVCENT192-12 | 12mm ø |
| PVCENT192-14 | 14mm ø |
| PVCENT192-16 | 16mm ø |
| PVCENT192-18 | 18mm ø |

| Code | Description |
|--------------|-------------|
| PVCENT192-20 | 20mm ø |
| PVCENT192-22 | 22mm ø |
| PVCENT192-24 | 24mm ø |
| PVCENT192-26 | 26mm ø |
| PVCENT192-28 | 28mm ø |
| PVCENT192-30 | 30mm ø |

Mitchell Trimmer



| Code | Description |
|------------|-------------|
| SC-ENT/013 | 165mm |

Miscellaneous

Lake Seeker

Double ended, for vestibule



| Code | Description |
|------------|-------------|
| SC-ENT/040 | 184mm long |


Negus Pusher

"Knot tier", metal handle




| Code | Description |
|------------|-------------|
| SC-ENT/165 | 203mm long |

Lempert Rugine



| Code | Description |
|-------------|--------------------|
| SC-ENT/064a | Narrow, 156mm long |
| SC-ENT/065a | Broad, 171mm long |

Stanford Cade Insu ator



| Code | Description |
|------------|-------------|
| SC-ENT/188 | 175mm long |

Irwin-Moore Tonsil Needle

203mm long



| Code | Description |
|-------------|-----------------|
| SC-ENT/164a | Curved to right |
| SC-ENT/164a | Curved to left |

Tubing Mount

For use with catheters, Luer



| Code | Description |
|------------|-------------|
| SC-ENT/016 | 25mm |